




Instructions for Use

CHEETAH® XTREME &
CHEETAH® XTEND

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	11
ES	Instrucciones para el uso	14
IT	Istruzioni per l'uso	17
DA	Brugsanvisning	20
SV	Bruksanvisning	23
NL	Gebruiksaanwijzing	26
PT	Instruções de Utilização	29
PL	Instrukcja użytkowania	32
JA	取扱説明書	35

5

6

CHEETAH® XTREME & CHEETAH® XTEND

Cheetah System

- Custom foot model for sprinting and recreational running.
- Maximum build height
 - Cheetah Xtreme 530mm (20 14/16")
 - Cheetah Xtend 451mm (17 3/4")
- Minimum build height
 - Cheetah Xtreme 375mm (14 3/4")
 - Cheetah Xtend 330mm (13")
- Maximum weight 147 kg (325lbs)
- One foot size

LAMINATION CONNECTOR (Figure 1)

FSX50009 (1-9)

- Lamination connector is laminated directly into the posterior socket wall.
- Placement of the lamination connector is recommended in the middle 1/3 of the socket for easier height adjustment.
- Make sure the socket surface is flat over the lamination connector to ensure good load distribution on the foot module when connected.

Bolts go through blade into lamination connector in socket. Apply Loctite 410 to all threads and torque to following specifications:

Category	Bolt Size	Torque
1-9	3/8"	45 Nm (33 ft-lbs)

The Cheetah® Xtreme and Cheetah® Xtend may also be connected by direct lamination.

FRICTION PADS

Always use friction pads between foot module and socket.

ALIGNMENT ADAPTER FSX50007 (Figure 2&3)

IMPORTANT: Only use alignment adapter for dynamic alignment adjustments.

- Determine adapter to socket position according to athlete height requirements.
- Attach alignment adapter to check socket using epoxy or bonding compound.
- Reinforce attachment using Tech Form Cast tape.
- Use adjustment screws to adjust foot contact point following alignment recommendations:

– Figure 2

1. Inversion/eversion foot adjustment
2. Height adjustment

– Figure 3

1. Dorsi/plantar flexion adjustment screw
2. Internal/external rotation adjustment screw

BENCH ALIGNMENT

Sagittal plane (Figure 4)

- Determine vertical socket position in alignment jig.
- Load line bisects socket and falls through the ground contact point or approximately 63mm (2.5") posterior to end of toe.
- Clamp foot into a 7° flexion position.

Coronal plane

- Introduce appropriate socket abduction/adduction using a jig.
- Determine appropriate internal or external socket rotation.
- Outset the foot module approximately 12 mm (1/2").
- Achieve proper height approximately 38mm (1 1/2") taller than the sound side to compensate for vertical compression.

NOTE: The recommendations for bench alignment are used as a starting point. Definitive alignment will be achieved during dynamic assessment using the alignment adapter.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

Cheetah Xtreme & Cheetah Xtend

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
WEIGHT LBS	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-325
Extreme Impact Level	1	2	3	4	5	6	7	8	9

CUTTING FOOT MODULE

- Ensure athlete has enough clearance during knee flexion.
- Trim foot module to the top end of the laminating connector.
- Cut foot module using a hack saw or band saw.

NIKE SPIKE PAD (Figure 5 & 6)

- Sand the bottom surface of the foot module
- Bond the Nike Spike Pad using contact cement or industrial strength glue
- Follow glue manufacturer bonding instructions

WATERPROOF

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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