
Instructions for Use

CHOPART

CHP00xyy


		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
DA	Brugsanvisning	19
SV	Bruksanvisning	22
NL	Gebruiksaanwijzing	25
PT	Instruções de Utilização	28
JA	取扱説明書	31

1


6 mm (1/4")
10 mm (3/8")


2


3


4

5

6

CHOPART FOOT PLATE (Figure 1)

- Uses a single carbon plate for direct lamination
- Build height 6 mm (1/4")
- Available in sizes 23–30 cm
- Not available in split toe version

BENCH ALIGNMENT (Figure 2 & 3)

- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- The load line taken from the center of the socket should fall through the midline of the foot module in a neutral M-L position. Depending on proximal or distal bearing socket design, use appropriate load line.
- Set appropriate toe rotation.
- Divide the foot into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions.

CHOPART INITIAL FITTING PROCEDURE

Preparation

From bench alignment, chopart foot plate should be bonded to the distal socket with potting epoxy or Siegelharz.

- Roughen carbon fiber surfaces (upper plate and socket) and clean with acetone to ensure best adhesion of epoxy or lamination.
- While in alignment jig, raise socket on the plate and lower unit back in place.
- Fill in around the outside edges with potting epoxy to create a smooth transition from the socket to the plate.
- Wrap with synthetic casting tape.
- Carry out dynamic alignment and make changes as necessary.
- During dynamic alignment you may cut off the tip of the heel and replace with soft foam material (**Figure 4**).

Dynamic Alignment

The roll over action of the foot can be influenced by the heel and the shape of the sole.

- Secure achieved dynamic alignment using a jig.
- Remove sole materials and synthetic casting tape and proceed to direct lamination.

Lamination Procedure

- The lay-up sequence involves posterior and M-L positioned carbon fiber cloth.
- Apply on posterior layer, then follow with a M-L layer.
- Continue to alternate using the required number of layers. The long posterior layer is folded in half over the M-L layer and reflected backup to the socket to form a T layup, locking the layers together and producing a clean lamination edge (**Figure 5**).
- Seal the not laminated part of the plates with tape.
- Pull a PVA bag over the prostheses and laminate with acrylic or epoxy resin.
- After lamination has set, sand off any protruding edges.
- Do not cut into the lamination.

- The standard foot cover may be modified to fit the prosthesis.
- For cosmetic finish the prosthesis can be covered with a 2 mm Pelite, waterproof sealed and color coated.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Low Impact Level	1	1	2	3	4	5	6	7	8
Moderate Impact Level	1	2	3	4	5	6	7	8	Special Order
High Impact Level	2	3	4	5	6	7	8	Special Order	Special Order

SUGGESTED REINFORCEMENT (Figure 6):

		Minimum* Layers by Category		
Layup Chart Sequence	Cloth Size	1-3	4-6	7-9
M-L Layers	76 mm × 355 mm 3" × 14"	3	4	6
Posterior Layers	76 mm × 355 mm 3" × 14"	2	2	3
Circumferencial Wrap	114 mm 4½"	2	2	3

*The appropriate number of layers depends on user activity and weight and must be evaluated by the CPO.

WATERPROOF

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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