
Instructions for Use







---

**FLEX-FOOT® JUNIOR**

Product number: JRPoxyyz

|    |                           |    |
|----|---------------------------|----|
|    |                           | 3  |
| EN | Instructions for Use      | 4  |
| DE | Gebrauchsanweisung        | 7  |
| FR | Notice d'utilisation      | 11 |
| ES | Instrucciones para el uso | 14 |
| IT | Istruzioni per l'uso      | 17 |
| DA | Brugsanvisning            | 20 |
| SV | Bruksanvisning            | 23 |
| NL | Gebruiksaanwijzing        | 26 |
| PT | Instruções de Utilização  | 29 |
| PL | Instrukcja użytkowania    | 32 |

# ENGLISH

---

## **FLEX-FOOT® JUNIOR (Figure 1)**

- Low profile foot module
- Build height: 46mm (1 13/16")
- Pyramid adapter 45-50 mm
- Maximum weight 45 kg (100 lbs)
- Available in 10 mm (3/8") heel height
- Available in sizes 16-21 cm

**Note:** Select foot category according to weight.

## **BENCH ALIGNMENT (Figure 2+3)**

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram (the mark on the brim of the foot cover matches the one third portion).

## **DYNAMIC ALIGNMENT**

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

## **CATEGORY SELECTION CHART**

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

| Weight Kg                | 15-24 | 25-34 | 35-45 |
|--------------------------|-------|-------|-------|
| Weight LBS               | 33-53 | 54-75 | 76-99 |
| Low to High impact level | 1     | 2     | 3     |

## **HEEL STIFFNESS**

### **Heel too soft**

#### *Symptoms*

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

#### *Action*

- Add polyurethane heel wedge.
- Shift socket anterior (or foot posterior).

**Heel too hard**

*Symptoms*

- Rapid heel to toe movement
- At initial contact the amputee has poor control of his prosthesis
- Minimal energy return feeling
- Knee may become unstable

*Action*

- Replace heel wedge with a smaller one
- Shift socket posterior (or foot anterior).

**HEEL WEDGES CM150011 (Figure 4)**

The wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The 2°, 4° & 6° wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics. A combination of wedges may be used

1. 2° – yellow
2. 4° – red
3. 6° – blue

**Temporary Wedge Placement**

- Cut the wedge to the width of the foot module.
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module.

**Note:** Do not force the wedge into the junction of heel and foot module.

- Secure in position bond to upper surface (**Figure 4**) with tape wrapped around the foot module (**Figure 5**).

**Permanent Wedge Placement**

- Apply adhesive on the upper side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

**HEEL ASSEMBLY**

If it is necessary to replace the heel component or heel bolts use Loctite 410 and torque to following specifications:

| Category | Bolt Size | Heel Bolt Torque |
|----------|-----------|------------------|
| 1-3      | M6        | 14 Nm (10ft-lbs) |

**Note:** A biannual inspection of attachment bolts is recommended: Check for corrosion and sign of wear. Replace bolts if needed.

**SPECTRA SOCK FCX63008**

To protect the cover and avoid noise a spectra sock is fitted over the foot module before installing the cover. The sock should be loose in the heel area and secured proximally with the cable tie provided.

**FOOT COVER JRX50YYZ**

Ensure the cover is completely inserted and the heel clip locks the foot module into the cover. Seal the cover to avoid debris entering.

**SHOE HORN CM100068**

When removing and installing the cover use shoehorn to avoid damage to foot or cover **(Figure 6)**.

**WEATHERPROOF**

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

**LIABILITY**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

**COMPLIANCE**

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee’s activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg ☆)

☆) Body mass limit not to be exceeded!  
For specific conditions and limitations of use see manufacturer’s written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases,

the test loads have been adapted adequately on the basis of the specified load level.

| ISO10328 |             |             |    |    |    |
|----------|-------------|-------------|----|----|----|
| Category | Weight (Kg) |             |    |    |    |
| 1        | 24          | ISO 10328 - | P3 | 24 | kg |
| 2        | 34          | ISO 10328 - | P3 | 34 | kg |
| 3        | 45          | ISO 10328 - | P3 | 45 | kg |

- This manual is intended for use by a certified prosthetist.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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