
Instructions for Use

FLEX-RUN™


with Nike Sole

Product number: FSX008xyy


		3
EN	INSTRUCTIONS FOR USE	4
DE	GEBRAUCHSANWEISUNG	7
FR	NOTICE D'UTILISATION	10
ES	INSTRUCCIONES PARA EL USO	13
IT	ISTRUZIONI PER L'USO	16
DA	BRUGSANVISNING	19
SV	BRUKSANVISNING	22
NL	GEBRUIKSAANWIJZING	25
PT	INSTRUÇÕES DE UTILIZAÇÃO	28
JA	取扱説明書	31

1


2


BENCH ALIGNMENT


3


4

5

6

ENGLISH

FLEX-RUN™ SYSTEM (Figure 1)

- Custom foot for recreational sports activities
- Preferred long distance running foot
- Build Height: 261mm (10 1/4"), clearance: 208mm (8 3/16")
- Maximum weight 130 kg (285 lbs)
- Endoskeletal attachment options

CONNECTOR OPTIONS

- Pyramid Adapter FWX00030 (1–6)
- Female Pyramid CM 150161 (1–6)
- Female Pyramid CM 150162 (7–9)
- Attachment plate FWX40001 must be added when using male or female adapter.
- Roughen top surface of Flex-Run (**Figure 2**) and attachment plate FWX40001
- Mix Siegelharz and apply to both surfaces
- Bond both parts together and clamp securely using a greased center bolt and the connector part
- Disassemble after resin has set and trim off excess
- Clean parts, assemble and torque using Loctite 410 to following specifications:

Center bolt	M10	41 Nm (30 lbs)
	M12	43 Nm (32 lbs)
Pyramid angulation	M 8	15 Nm (9 lbs)

BENCH ALIGNMENT (Figure 3 & 4)

SAGITTAL PLANE

- Place the ground contact point of the foot module 25–40 mm (1–1½") in line with load line.
- Determine the appropriate socket flexion according to amputee requirements.

CORONAL PLANE

- Introduce appropriate socket abduction/adduction.
- Allow for up to 5° of external rotation.
- Achieve proper height approximately 25–50 mm (1–2") taller than the sound side with shoe to compensate for vertical compression.
Note: The recommendations for bench alignment are used as a starting point. Definitive alignment will be achieved during dynamic assessment.

DYNAMIC ALIGNMENT

Physical ability and length of residual limb will affect alignment. Best performance will be achieved with best achievable symmetry of gait. Each athlete will have her/his own optimal alignment, which depends on the individual's running style. Optimal alignment is critical in order to reduce energy expenditure and eliminate excessive muscle and ligament strain while running:

- Consult with your Össur clinician before using the running foot for the first time.
- Check for best symmetry of gait.
- Stride length can be affected by A-P alignment and vertical deflection socket to foot.
- Vertical deflection is between 25 mm (1") and 50 mm (2") on average.

SOLE ATTACHMENT (Figure 5 & 6)

- Slide the Nike sole on to the foot and make sure all the clips go over the foot.
- Pull the rubber strap of the sole over the medallion and slide it under. When detaching the sole, start by pulling the strap away from the medallion, then slide the sole off the foot.
- For tighter fit of the sole or to reduce noise bond insert inside sole using Loctite, keep under pressure for 10 minutes cure

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	37-44	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130
WEIGHT LBS	81-96	97-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287
High to Extreme Impact Level	1	2	3	4	5	6	7	8	9

TROUBLESHOOTING

Make sure socket fit and suspension are ideal. Symmetry in running should be assessed in the optimal individual speed phase.

If the duration of loading the prosthesis is longer than the contralateral side:

- Foot module is too soft
- Prosthesis is too short
- Foot is too far anteriorly aligned

Duration of loading the prosthesis is shorter:

- Foot module is too stiff
- Prosthesis is too long
- Foot is too far posterior aligned

TRANSFEMORAL ATHLETE CONSIDERATIONS

Moving the foot in A-P position will affect hip flexion, stride length and trunk motion.

WATERPROOF

A waterproof device allows use in a wet and/or humid environment and temporary (<30 mins.) submersion of up to 1m in fresh water. Exposure to fresh water under these circumstances shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

The components used with the foot need to be waterproof as well.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is

not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass.

In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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