
Instructions for Use


VARI-FLEX® JUNIOR


		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
EL	Οδηγίες χρήσης	37
NL	Gebruiksaanwijzing	42
PT	Instruções de Utilização	46
PL	Instrukcja użytkowania	50
CZ	Návod k použití	54
TR	Kullanım Talimatları	58
RU	Инструкция по использованию	62
JA	取扱説明書	67
ZH	中文说明书	71
KO	사용 설명서	75

1


Size	Build height
18	120 mm (4 3/4")
21	123 mm (4 13/16")
24	124 mm (4 7/8")


2


3


4


5

6

7

VARI-FLEX® JUNIOR FOOT SYSTEM (Figure 1)

- Minimal Clearance 120mm-124mm (4 3/4" - 4 7/8")
- Maximum weight 55 kg (121 lbs)
- Available in sizes 16-24 cm

SAFETY PRECAUTIONS

WARNING: Statement that alerts the user to the possibility of serious injury, or serious adverse reactions associated with the use or misuse of the device.

CAUTION: Statement that alerts the user to the possibility of a problem with the device associated with its use or misuse.

PRECAUTION: Please follow the safety precautions as failure to do so can lead to malfunction of the device and risk injury to the user. Prosthetists should make sure their clients are aware of, and understand, the safety precautions.

WARNING: Make sure that the user knows if he/she feels or notices any change of function or alignment of the device, he should contact his healthcare practitioner. This noticeable function loss may be the result of reduced stiffness or loss of support.

CAUTION: Make sure that the user knows that he needs to replace socks that have holes in it because otherwise it may damage the foot cover.

PRECAUTION: The foot is designed to be used with spectra socks and foot cover to prevent noise and for proper shoe fit.

PRECAUTION: The device is weather-proof but not resistant to corrosion. Therefore, the device should not come into contact with salt water or chlorinated water. The device is not designed to be used under extreme conditions like diving or jumping into water. If for some reason the foot components get wet parts should be dried with a lint free cloth. Should the device come into contact with salt water, chlorinated water, dust or dirt it must be rinsed with fresh water immediately.

INDICATIONS FOR USE

The foot is a non-invasive, single patient, reusable prosthetic device. The foot is intended for unilateral and bilateral lower extremity amputees.

CONTRAINDICATIONS FOR USE

N/A

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	15-24	25-34	35-45	46-55
WEIGHT LBS	33-53	54-75	76-99	100-123
16	1	N/A	N/A	N/A
17	1	N/A	N/A	N/A
18	1	2	N/A	N/A
19	1	2	N/A	N/A
20	N/A	2	3	N/A
21	N/A	2	3	N/A
22	N/A	2	3	4
23	N/A	2	3	4
24	N/A	N/A	N/A	4

BENCH ALIGNMENT (Figure 2)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in figure 2
- The 1/3 position is marked on the inside wall of the foot cover

DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by:

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff
- Knee may hyperextend.

Actions

- Add polyurethane heel wedge
- Shift socket anterior (anterior or foot posterior).

Heel too hard

Symptoms

- Rapid heel to toe movement.
- At initial contact the amputee has poor control of his prosthesis.
- Minimal energy return feeling.
- Knee may become unstable.

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HEEL WEDGES (Figure 3)

A wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge.

Temporary Wedge Placement

- Cut the wedge to the width of the foot module
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module
- Secure in position with tape wrapped around the foot module
(Figure 4).

Permanent Wedge Placement

- Apply adhesive on the lower side of the wedge only (1).

- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module (**Figure 5**).

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds.

HEEL ASSEMBLY

If it is necessary to replace the heel component use Loctite 410 and torque to following specifications:

Category	Bolt Size	Heel Bolt Torque
1-4	1/4"	14Nm (10ft-lbs)

HEEL DIVIDER

The heel divider should be placed in the last third of the splitted heel part. To secure its position a drop of instant adhesive can be used on one side.

BOLT CAP REMOVAL

If it is necessary to remove the bolt cap locate the foot module/bolt cap junction and use a sharp knife to separate.

FLEX-FOOT SOCK

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area and secured proximally with the cable tie provided.

SHOE HORN (Figure 6 & 7)

When removing and installing the cover use shoehorn to avoid damage to foot or cover.

MAINTENANCE

The device is designed for low-maintenance. The device should be inspected every six months for signs of unusual wear, by a prosthetist familiar with the product.

WEATHERPROOF

A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

USAGE ENVIROMENT

- Operating temperature: -15°C (+5°F) to 50°C (122°F).
- Operating humidity: 0%- 100%. Relative humidity.
- Shipping and storage temperature: -15°C (+5°F) to 50°C (122°F).
- Storage humidity: 0%- 100%. Relative humidity (non-condensing).

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee’s activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer’s written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been

adapted adequately on the basis of the specified load level.

ISO10328					
Category	Weight (Kg)	Lable text			
1	24	ISO 10328 -	P3	24	kg
2	34	ISO 10328 -	P3	34	kg
3	45	ISO 10328 -	P3	45	kg
4	55	ISO 10328 -	P3	55	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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