






Instructions for Use






TOTAL KNEE® JUNIOR

Product Number: 1100

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	11
FR	Notice d'utilisation	18
ES	Instrucciones para el uso	25
IT	Istruzioni per l'uso	32
DA	Brugsanvisning	39
SV	Bruksanvisning	45
EL	Οδηγίες χρήσης	51
NL	Gebruiksaanwijzing	58
PT	Instruções de Utilização	65
PL	Instrukcja użytkowania	72
JA	取扱説明書	79
ZH	中文说明书	85

ENGLISH

NUMBERING IN FIGURES

1. Alignment reference line
2. Socket
3. Knee centre reference
4. Front link
5. Balance unit
6. Pivot axis (alignment axis)
7. Housing
8. Back link
9. Bottom link
10. Chassis
11. Bumper
12. Shim
13. Extension promoter
14. Spring cup
15. Spring
16. Ferrule
17. Cable
18. Retaining Ring
19. Spring cup seat
20. Knee cap – small
21. Knee cap – large
22. Shin ferrule
23. Bearings (lubrication points)
24. Release points

FUNCTIONS AND BENEFITS

- Increased knee stability
- Increased toe clearance in swing phase
- Lightweight design
- 160° knee flexion – suitable for sitting, kneeling, cycling and other activities where greater knee flexion is desired
- 24-month warranty

INDICATIONS FOR USE

- Junior amputee up to 45kg (100lbs) patient weight
- The Total Knee is exclusively used for lower limb prosthetic fittings

CLEARANCE (Figure 1)

ALIGNMENT AND ADJUSTMENT (Figure 2)

ADAPTER OPTIONS (Figure 3)

- 4- prong socket adapter (AJ-114040).
- Pediatric 4-hole adapter (AJ-122300) (1155).
- Pediatric pyramid adapter (AJ-834300) (1125).

Important: Always use original Total Knee/Össur adapters. Fully engage threads. Align pinch bolt 90° anteriorly for male pyramid, prong adapters and Euro 4-hole adapter (Figure 3). Threads are cut for maximum strength in this position.

Adapters from other manufacturer may not fully engage all threads. They may compromise strength and void the warranty. Never tighten adapters against the head.

When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the knee.

Important: Contact between socket/adapter and proximal 2/3 of back link must be avoided when fully flexing knee (**Figure 4**). Contact in this area can damage knee and void the warranty. Contact between socket and back link below balance unit is acceptable.

If pylon is too short for changed setup, do not compensate by putting a spacer into the knee tube adapter as this may compromise strength and void the warranty (**Figure 5**). Use a new Össur pylon cut to correct length for maximum patient safety.

BENCH ALIGNMENT (Figure 2)

- Assemble prosthetic foot, pylon and knee.
- Establish knee centre height, measured from the knee center reference(3). Tighten tube clamp to 8 Nm(6 ft-lbs) (**Figure 6**)
- Place appropriate heel height under heel or fit patient's shoe.
- Establish vertical pylon for proper knee function. Alignment reference line (1) from pivot axis (6) should pass through posterior 1/3 mark on Flex-Foot® cover (or see other foot system instructions).
- Attach socket with selected adapter combination to knee. Tighten to torque as specified in adapter IFU. Ensure flexion/extension and abduction/adduction angles are as determined in patient assessment. Alignment reference line (1) from bisection of socket on ischial level (2) should pass through pivot axis (6) of knee and posterior 1/3 mark on Flex-Foot cover.
- Make sure total length of prosthesis is correct.

NOTE: If alignment reference line (1) is too far anterior to pivot axis(6), more effort may be required to initiate knee flexion.

BUMPER SELECTION (Figure 7)

The bumper determines the amount of stance flexion.

- Insert small flathead screwdriver under bumper and lift upwards.
- Remove bumper (leave the shims under it in place).

Patient Weight	Colour	Type
14-20 kg (30-45 lbs)	Green	Soft
20-29 kg (45-65 lbs)	Black	Regular
29-36 kg (65-80 lbs)	Red	Firm
36-45 kg (80-100 lbs)	Yellow	Extra Firm

- Install appropriate bumper according to patient weight.
- Active patients may require a stiffer bumper
- Patients who have previously worn a prosthesis without a stance-flexion feature may initially prefer a stiffer bumper.
- We suggest trying a softer bumper after two weeks of trial. This helps the patient develop security with the knee. Doing this will also help make the transition to the stance-flexion characteristic more acceptable.

STATIC ALIGNMENT

For safety, please make initial adjustments with patient standing between parallel bars!

- Fit prosthesis and check for correct length.
- Make sure foot is flat on the floor, knee is in neutral position and socket angles are correct.
- When the patient shifts weight onto the prosthesis, the geometric lock is activated and the bumper should be slightly compressed. Modify alignment if necessary.
- Explain function of stance flexion/geometric lock to patient. Ask patient to activate it by stepping forward with prosthesis and transferring weight onto heel.
- Patient should become accustomed to activating geometric lock and resulting stance flexion movement.
- Knee cannot collapse as long as geometric lock is activated.

DYNAMIC ALIGNMENT

Swing Control Adjustments

The pediatric Total Knee is designed to operate with free extension resistance to provide the quick response required by most children. Therefore, swing phase adjustments are completed through use of the extension promoter and extension assist.

Bumper adjustment

- Try different bumpers and select the most comfortable one.
- Final bumper selection may vary from bumper indicated for patient's weight according to table.
- Soft bumpers increase the locking moment of the Geometric Lock and increase stance flexion motion.
- Hard bumpers decrease the locking moment of the Geometric Lock and decrease stance flexion motion.
- Foot selection influences bumper selection.
 - a. For soft heels, choose a softer bumper.
 - b. For hard heels, choose a hard bumper.

Shim adjustment (Figure 7)

Adding shims changes the release point of the Geometric Lock allowing for swing phase knee flexion. This affects knee stability and ease of walking and sitting. The shims are thin plastic plates which are installed underneath the bumper to adjust its height. A change in height affects the knee's release point.

- More shims: easier release into swing flexion and decreased stance flexion
- Less or no shims: increased safety, increased stance flexion and delayed release into swing flexion

NOTE: There are two different types of shim. The thicker ones are 0.76mm thick and clear, the thinner ones are 0.38mm thick and grey. We do not recommend using more shims than those two combined. Adding more shims can eliminate the Geometric Lock and affect stability of the knee during stance.

Gradually increase shim thickness until the user feels premature release of the knee. Reduce shims until the user finds a balance between security and release point. Please make initial adjustments inside parallel bars.

NOTE: If no shims are used, the knee will exhibit some flexion/extension play during stance. If the patient finds this motion unacceptable, this can be reduced by adding shims.

EXTENSION PROMOTER (Figure 7)

The extension promoter (13) reduces excessive heel rise and promotes faster knee extension. It is set at minimum when shipped from Össur.

- Flex the knee to access the extension promoter adjustment screw.
- To reduce excessive heel rise, especially when walking fast, turn the extension promoter adjustment screw clockwise with a large flathead screwdriver.

EXTENSION ASSIST (Figure 9)

The extension assist is an **OPTIONAL** feature that helps get the knee back to full extension. We recommend use of the extension assist only when the amputee is unable to extend the knee to engage the geometric lock.

- Beneficial for users having difficulty achieving full extension.
- Can help during fast walking.
- Brings knee to full extension consistently.
- Places positive locking bias on Geometric Lock to enhance stability.

Available with two different springs (yellow=regular, blue=firm).

INSTALLATION (Figure 9)

- Remove bumper and shims (**Figure 7**)
- Install spring cup (14) into hole inside 30mm tube clamp of knee (19)
- Firmly seat spring cup by gently tapping with rubber mallet
- Place ferrule (16) over cable (17)
- Place spring (15) over cable (17) and ferrule (16) (start with the yellow spring)
- Insert cable (17) and spring (15) into spring cup (14), ball end of cable first
- Slide ball end of cable through hole in knee (19) and insert ball into the slot on bottom link (9) using small screwdriver or needle nose pliers.
- Reinstall the selected shims (12) and bumper (11).

FOLLOW-UP CHECK

- Reassess adjustments after a couple of weeks
- It may be necessary to change shim and bumper adjustments after user has gained some experience with the knee
- Increased experience, security and voluntary control may require change to a harder bumper and/or increased shim thickness.

COSMETIC OPTIONS (Figure 8)

We recommend a discontinuous cosmetic solution which does not interfere with knee function. Continuous cosmetic solutions are also possible.

Discontinuous Cosmetics

Ref. No	Description	Notes
1145	Junior Shin ferrule	Protects unit from foam cover - maintains space for movement within foam cover
1150	Junior Knee cap - small	Use to prevent tearing at the knee
1151	Junior Knee cap - large	Use for patients who kneel often

MAINTENANCE

The Total Knee 1100 is designed for low-maintenance. It is warranted against defective materials and workmanship for 24-months. We recommend:

- Total Knee 1100 be inspected every six months for signs of unusual wear, by a prosthetist familiar with this product.
- Stance flexion bumper must be replaced if cracked or worn, depending on patient's activity level.
- If knee is subjected to excessive moisture or corrosive environment, it's recommended to clean and lubricate the knee frequently.

CLEANING

- Wipe knee with soft cloth moistened with small amount of kerosene. DO NOT use solvent stronger than kerosene. If kerosene is not available, use general purpose oil (3in1) or sewing machine oil.
- DO NOT dip knee or pour solvent over knee. Bearings and seals may be damaged.
- DO NOT use compressed air to clean knee. Air forces pollutants into bearings and may cause malfunctions and wear.

LUBRICATION (Figure 10)

Three sets of roller bearings are visible. Place a few drops of machine oil on the roller bearings (general purpose oil (3in1) or sewing machine oil can also be used). Move the knee several times and wipe off excessive oil with a soft cloth. No other parts need external lubrication.

PRODUCT-SPECIFIC TRAINING

Instructing new users is essential to achieve a successful rehabilitation. For safety, initial adjustments and gait training should be done with patient standing between parallel bars.

LEVEL GROUND WALKING

- Explain to the user the function of stance flexion/geometric lock. Patients who previously wore a knee without stance flexion function may fear collapse of knee during loading response. They must be trained to allow the knee to move into stance flexion and not to prevent it by strong hip extension
- Explain the stance release mechanism to the patient. Patient's centre of mass must be anterior to the pivot axis of knee. In this position, very little hip flexor effort is required to initiate knee flexion. Knee flexion will be initiated through rotation of the hip as weight is transferred onto the ball of the foot during terminal stance.
- Allow patient to stand in the parallel bars with their weight over the ball of the foot. Initiate knee flexion with pelvic rotation and slight hip flexion. Repeat several times.

- Start walking within parallel bars. As confidence increases, start walking outside of the parallel bars.

SITTING DOWN

- Have the patient move their centre of mass forward over the toe
- Place feet evenly (or prosthetic foot slightly posterior), shift weight anterior onto forefoot. Knee will release
- Use sound limb to accept weight for sitting

STAIR DESCENT

Walking step over step down stairs with Total Knee requires much practice and should only be attempted by trained Total Knee users. Begin practicing on lowest step and use handrail.

- Place entire foot on step.
- Move centre of mass over forefoot as for sitting. Knee flexion will be initiated. Knee flexion will be rapid. Use caution!
- Increase number of steps as patient becomes confident.

Not all users will be able to learn step over step stair descent.

MANUAL FLEXION

When sitting with fully extended Total Knee it may be necessary to flex knee manually. Press two points shown in **Figure 11** to release knee into flexion. Be careful not to pinch fingers in the linkages. (Patient may also reach out and pull or tap toe of prosthesis to release.)


CAUTION! Avoid placing hands or fingers near moving joints.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to three million load cycles. Depending on the amputee's activity this corresponds to a duration of use of three to five years. We recommend carrying out regular yearly safety checks

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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