
Instructions for Use


FLEX-SYMES™

Product Number : SLPOXYYO


		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	11
FR	Notice d'utilisation	17
ES	Instrucciones para el uso	23
IT	Istruzioni per l'uso	29
DA	Brugsanvisning	35
SV	Bruksanvisning	41
NL	Gebruiksaanwijzing	47
PT	Instruções de Utilização	53
JA	取扱説明書	59

1


Size **Build height**
 22-30 cm = 60 mm (2 5/16")








BENCH ALIGNMENT

2

3

FLEX-SYMES™ COMPONENTS (Figure 1)

- Circular alignment wedges (SLX3001X) allowing up to 6° angulation and adjustments in AP/ML planes
- Available in sizes 22-30cm.
- Minimal clearance:
 - 53 mm for sizes 22-24
 - 55 mm for sizes 25-27
 - 57 mm for sizes 28-30

Figure 1

1. Foot Module
2. Heel (SLX2XYY0)
3. Heel fastener hardware kit (CM150004)
4. Socket Adapter (SLX5000X)
5. Alignment Disc Kit (SLP30001)
6. Central Bolt (SLX3002X)
7. Anti-Rotation Bolt Kit (SLX3000X)
8. Heel Wedge Kit (CM150011)

BENCH ALIGNMENT (Figure 2 & 3)

- Attach heels to the foot module without glue, fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction /adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram (the mark on the brim of the foot cover matches the one third portion).

FABRICATION PROCEDURES (Figure 4)

- Leave foot assembled to the socket adapter and protect foot with appropriate film/sheet barrier.
- Roughen attachment surface, position socket in the desired alignment and fill the space beneath the socket with epoxy putty (CM160009).
- Let it harden (3-5 minutes)
- Wrap synthetic casting tape and proceed with dynamic alignment for temporary use only.

Note: Do not use polyurethane foam build-up since this compresses and can lead to bolt or lamination failure.

IMPORTANT

Fasten the attachment bolts through the heel. Make sure there is no interference with the central bolt head when heel is compressed. Verify torque specifications. Position heel divider (CM140003) as demonstrated in picture (**Figure 8**).

DYNAMIC ALIGNMENT

Alignment Wedges

Use the alignment options of the circular alignment wedges to fine tune the dynamic alignment. The alignment wedges provide angular adjustments up to 6 degrees.

Instructions:

- Two 3mm steel pins are provided in the kit.
- Insert pins into alignment wedges.
- Twist alignment wedges against each other to fine tune alignment.

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-356
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	N/A
High Impact Level	2	3	4	5	6	7	8	9	N/A	N/A

HEEL FUNCTION AND ADJUSTMENTS

The heel stores energy after initial contact, thus slowly releasing it at midstance.

The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by:

- A-P positioning of foot.
- Dorsi-Plantarflexion.
- Heel Stiffness.

Heel too Soft

Symptoms

- Foot will come to foot flat position too early, (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

Action

- Add polyurethane heel wedges.
- Shift socket anterior/ foot posterior.

Heel too Hard

Symptoms

- Rapid heel to toe movement.
- Amputee has poor control of prosthesis.
- Toe will feel too stiff, causing minimal energy return.
- Knee may become unstable.

Action

- Replace heel wedge with smaller wedge.
- Shift socket posterior or foot anterior.

HEEL WEDGES CM150011 (Figure 5)

The wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The 2°, 4° & 6° wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics. A combination of wedges may be used.

1. 2° – yellow
2. 4° – red
3. 6° – blue

Temporary Wedge Placement

- Cut the wedge to the width of the foot module.
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the foot/heel junction.
- Secure in position, bond the wedge to upper surface with tape wrapped around the foot module.

Permanent Wedge Placement

- Apply adhesive on the upper side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module.

LAMINATION TECHNIQUE

- Mark alignment position after dynamic alignment and remove foot module from socket.
- ***Note:** If alignment wedges are straight, both wedges may be removed to decrease clearance and complexity! This will lower the height by 6mm.*
- Invert the socket and clamp securely in the lamination jig.
- Clean the groove of the aluminum socket adapter with acetone or similar to ensure maximum bonding of the fibers to the aluminum.
- Roughen attachment surfaces and seal center hole (by using plasticine or similar).
- Apply appropriate carbon lay-up (see table 1). Tie off the overhanging braids around the base grove using a heavy thread. Ensure that the tie stays tight.
- Finish with a stockinette for cosmetic finish, if desired.
- Pull a PVA bag and laminate with acrylic or epoxy resin.
- Tape with PVC-tape around the distal end.
- Tie off excess resin before it reaches the gel state. Allow resin to harden.
- Grind off the excess hardened resin around the edge of the aluminum socket adapter. Expose the attachment surface. Care must be taken not to grind the aluminum piece itself!

Table 1: Suggested Reinforcement

Sizes	Category	Layers of carbon fiber
22 – 24	1 – 6	Minimum 4*
25 – 30	1 – 6	Minimum 6*
25 – 30	7 – 9	Minimum 8*

*The appropriate number of layers depends on user activity and weight and must be evaluated by the CPO. It is recommended to use heavy carbon fiber braids.

ASSEMBLY (Figure 6, 7 & 8)

- Position circular alignment wedges according to dynamic alignment.
- Assemble and torque center bolt using Loctite 410, see table 2.

Table 2: Torque Specifications for central bolt.

Category	Torque
1 – 6	40 Nm (30ft-lbs)
7 – 9	60 Nm (44ft-lbs)

Note: Applies to all sizes.

- Guided by the positioning hole, drill a hole through the circular alignment wedge/s and the socket adapter. Tap threads for the anti-rotational bolt, see table 3 for specifications.

Table 3: Thread specifications for anti-rotation bolt.

Sizes	Drill ø [mm]	Tap
22 – 27	5.0	M6 x 1.0
28 – 30	6.8	M8 x 1.25

- Assemble, using Loctite 410, anti-rotation bolt (CM130326/7) and anti-rotation washer (SLX30055/6) to secure the alignment. Torque to 10 Nm, 8 ft-lb.
- Position the heel divider (CM140003) in the posterior slot of the splitted heel part, see fig. H. Bolt the heel to the foot module and torque to specifications, see table 4, using Loctite 410. Make sure there is no interference with the central bolt head when the heels are compressed. To secure the position of the heel divider a drop of instant adhesive can be used on one side.

Table 4: Torque specifications for heel attachment bolts.

Category	Bolt diameter [in]	Torque
1 – 6	1/4	14 Nm (10ft-lbs)
7 – 9	5/16	28 Nm (21ft-lbs)

Note: Applies to all sizes

IMPORTANT

It is important that the torque specifications are met.

COSMETIC COVER FCA0XXL/R

The foot section may be covered with the anatomical foot cover. Add "BR" to the part # for brown cover.

The Foot cover may need some trimming if the socket does not fit into the foot cover. A foot cover attachment cannot be used for FLEX-SYMES.

FLEX-FOOT SOCK FCX63006/7

To protect the cover and avoid noise, a Flex-Foot sock is fitted over the foot module before installing the cover. The sock should be loose in the heel area.

SHOEHORN CM100068 (Figure 9)

When removing and installing the cover, use shoehorn to avoid damage to foot or cover.

MAINTENANCE

The FLEX-SYMES can be used in water but warranty terms only apply if the following conditions are met:

- The foot is cleaned with tap water after use in salt water and / or freshwater.
- Immediately following contact with water, take the the foot cover off the foot and dry thoroughly.
- The components used with the foot need to be waterproof as well.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guideline.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee’s activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

ISO 10328 - “P” - “m”kg ☆)

⚠

*) Body mass limit not to be exceed!

⚠

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads

have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg
5	88	P4	ISO 10328 -	P5	88	kg

Category Össur high activity						
Category	Weight (Kg)	Lable text				
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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